

The logo features the word "JET" in a bold, white, sans-serif font, with a small "TM" trademark symbol to its upper right. Below "JET" is a smaller, white, lowercase "c" with a horizontal line passing through its center. This entire logo is contained within a teal circle. The background of the page is a white field with a pattern of green and teal circles of various sizes, arranged in a way that suggests a stylized, upward-pointing arrow or a series of overlapping paths.

JETTM
c

The word "cache" is written in a white, lowercase, sans-serif font. A small "TM" trademark symbol is positioned to the upper right of the letter "e". The text is centered within a solid teal horizontal bar that spans the width of the page.

cacheTM

READY TO USE JET? LET'S GET STARTED.

JET is a ClO₂ sporicidal foam for the disinfection of non-porous surfaces, delivering sporicidal, mycobactericidal, virucidal, fungicidal, yeasticidal, and bactericidal efficacy in 1 minute.

PRE-USE CHECKS

Some important reminders before starting: Use disinfectants (product type 2) safely. Always read the label and product information before use. Wear the appropriate personal protective equipment (PPE).

JET is for professional use only. Do not apply to brass, copper, aluminium, mild steel, or silver surfaces. Ensure the surface is free from heavy soiling before use.

To assemble JET, take a trigger, unscrew the caps on the bottles, and following the red A and blue B indicators on the bottle label, insert the tubes into the respective bottle necks. Fasten the levers anti-clockwise on each side of your trigger to secure. Ensure the nozzle is aligned to the top line on the trigger head. JET is now ready for use.

HOW TO USE



Apply at least 2 doses onto a wipe.



Wipe onto the surface.



Ensure the surface remains wet for 1 minute. Apply additional solution as needed. Leave to **dry**.

BIG IMPACT

When you can't recycle, reuse! Your JET trigger should be reused with the second JET bottle. To do this, unfasten the levers clockwise to remove the trigger. Before following the assembly instructions above, give the dip tubes a rinse with tap water and dry with a paper towel.

After use, empty JET bottles and caps can be sent for recycling.

HOW JET WORKS

JET incorporates two separate compartments that contain Base Solution (citric acid) and Activator Solution (sodium chlorite). When mixed by pulling the trigger, chlorine dioxide (ClO₂) chemistry is generated.

ClO₂ eliminates microorganisms by destroying their molecular structures through oxidation. This chemical process means that harmful microbes cannot build resistance to ClO₂.

STORAGE CONDITIONS

JET can be kept for 2 years before being used, and should be stored in a cool dry place, out of direct sunlight, between 10°C - 35°C.

Once opened, JET has an in-use life of 3 months. You can track this by noting the day you opened JET on the bottle label.



CONTENTS

- 400ml Activator Solution, sodium chlorite (2.5g/kg)
 - 400ml Base Solution, citric acid monohydrate (50g/kg)
 - Active at point of use: chlorine dioxide (ClO₂) (0.2g/kg)
- JET solutions are not classified as hazardous according to GB and EU CLP regulations.

This Biocidal Product Type 2 (Disinfectant) has a transitional status (no authorisation number) under the Biocidal Products Regulations EU BPR and GB BPR.

For supporting documents such as safety data sheets, microbiological test data and reports, contact Tristel, your local distributor or visit www.tristel.com.



TEMPERATURE LIMIT



KEEP AWAY FROM SUNLIGHT



DO NOT USE IF PACKAGE IS DAMAGED



CONSULT INSTRUCTIONS FOR USE



MANUFACTURER



USE BY DATE

LOT

BATCH CODE



cache™

Cache is a collection of products for the cleaning and sporicidal disinfection of surfaces brought to you by Tristel.



Tristel Solutions Ltd, Lynx Business Park, Fordham Road,
Snailwell, CB8 7NY, UK T +44 (0)1638 721500
E mail@tristel.com - W www.tristel.com

CCH/JET/050 - V1 - 01/OCT/2021

