

cache™

APPLICATORS

READY TO APPLY SOLUTION? LET'S GET STARTED.

These applicators are designed to apply Cache solutions onto non-porous surfaces and non-invasive medical devices. They can be paired with any of the 500ml bottles within the Cache collection.



The Doser delivers 5ml doses of liquid.



The Foamer delivers 1ml doses of foam.



The Sprayer delivers 1ml doses of liquid.

PRE-USE CHECKS

A few reminders before starting: These applicators are for professional use only and are not intended to be used with invasive medical devices.



All three applicators can be rinsed out and re-used. Be sure to replace your applicators after 1 year, or if they no longer dispense as expected.



Brought to you by:

Tristel Solutions Ltd, Lynx Business Park, Fordham Road, Snailwell, CB8 7NY, UK - T +44 (0)1638 721500
E mail@tristel.com W www.tristel.com



Tristel NV,
Smallandlaan
14 B, 2860
Antwerp,
Belgium

